


AMBULATORY INFUSION PUMPS

HOME PUMP C-SERIES

Soft shelled elastomeric pump with a fill port and tubing with an integrated particulate and air-eliminating filter.

There is a flow-restrictor. The most common sizes are :

Code	Nappi	Volume	Flow Rate	Duration
C060020	161807*001	60ml	2ml/hr	30 hours
C100020	161811*001	100ml	2ml/hr	2 days
C125050	161810*001	125ml	5ml/hr	1 day
C270010	161812*001	270ml	1ml/hr	11 days
C270020	161813*001	270ml	2ml/hr	5.5 days
C270050	161814*001	270ml	5ml/hr	2 days
C270100	161815*001	270ml	10ml/hr	1 day



The HOME PUMP ECLIPSE is available in volumes of 50—500ml at flow rates of 30 mins – 4 hours

MINI-PUMPS

Hard shelled pump with a return valved fill port and tubing with an 0.2 micron integrated particulate and air-eliminating filter. There is a flow-restrictor. The following sizes are available:

Code	Nappi	Volume	Flow Rate	Duration
DPI6020D	467012*001	60ml	2ml/hr	30 hours
DPI10020D	665696*001	100ml	2ml/hr	2 days
DPI10040D	665698*001	100ml	4ml/hr	1 day
DPI275100D	665699*001	275ml	10ml/hr	1 day
DPI27550D	665699*001	275ml	5ml/hr	2 days
DPI27520D	662700*001	275ml	2ml/hr	5 1/2 days



These pumps allow the patient to receive medication at a pre-determined rate without having to stay in hospital. They are often used for chemotherapy for oncology patients, anti-biotics for cystic fibrosis patients, for some enzyme disorders and for pain management.