


**DRESSINGS**

ZEBRA MEDICAL has a range of high quality dressings in different sizes, targeting the securement of a variety of different lines.

**SORBAVIEW (SV40xXT)**

Semi-permeable transparent polyurethane window dressing with an integrated absorbent border that wicks away any moisture. It is surrounded by a non-woven fabric. The skin is able to breathe, but is a barrier to liquids and bacteria. Can remain in position for 7 days. Boxes of 50

**VERSADERM (VSD30ST)**

Versaderm is a 6 x 8 cm transparent dressing with a window.

It allows the skin to breathe, but is impermeable to virus', bacteria and liquids.


Boxes of 100

**SECUREFIX PLUS (2 sizes)**

Sterile extra-thin adhesive polyurethane dressing .

5 x 7 cm (00022776000000) - boxes of 100

10 x 10 cm (00022786000000) - boxes of 50

**ZONIS SILVERSITE ANTIMICROBIAL ACCESS SITE DRESSING**

25mm diameter with 4mm hole (SS853XT)

37.5mm diameter with 7mm hole (SS171XT)

A slit circular shaped sterile dressing consisting of an absorbent polyurethane foam layer that provides for absorption of wound exudate and a second layer of ionic silver alginate matrix that provides broad anti-microbial effectiveness.

**EMO SOFFIX POST — BLOOD PLASTERS—38mm x38mm (022595000000)**

Small dressing ideal for use after blood, IV cannula or port needle removal.

The central pad absorbs 10 times it's weight, swelling up and making it an effective pressure dressing



ZEBRA MEDICAL is committed to supplying high quality dressings and “targeting zero” w.r.t. infections