

POLYPERF PORT NEEDLES—Straight, curved and winged

When accessing a port it is important to use a non-coring Huber needle to prevent damage to the self-sealing silicone of the port. Zebra Medical stocks a comprehensive range of needles — straight, curved, safety, positive pressure and CT.

STRAIGHT AND CURVED NEEDLES

These are ideal for bolus doses of medication or flushing a port.

- Straight Huber needle 22G x 30mm (512507)



POLYPERF WINGED HUBER NEEDLE WITH CONNECTING LINE (Boxes of 10)

- DEHP and Latex free
- Available with and without lateral injection site
- 90° Curved needle with luer lock and clamp
- Colour coded for easy gauge identification



Product code	Nappi codes	Without Y-site
582011	665727*001	G19 x 20 mm (3/4 “)
582511	665740*001	G19 x 25mm (1”)
583511	665735*001	G19 x 35mm (1,5”)
581709	665741*001	G20 x 17mm (1/2”)
582009	665730*001	G20 x 20mm (3/4”)
582509	665732*001	G20 x 25mm (1”)
583509	665743*001	G20 x 35mm (1,5”)
581707	665733*001	G22 x 17mm(1/2”)
582007	665739*001	G22 x 20mm (3/4”)
583507	665747*001	G22 x 35mm (1,5”)



Product code	Nappi codes	With Y-site
592011	462800*001	G19 x 20mm
592009	523106*001	G20 x 20mm
592509	523785*001	G20 x 25mm

The above are available in a variety of gauges and lengths.
Please contact us for any more information